

The device for enhancement  
of individual body defenses  
and correction of natural  
electromagnetic emission  
of natural biosystems

DeVita AP  
Base model

Methodological recommendations

## Introduction

You are holding a unique, unparalleled anywhere in the world, electromagnetic «DeVita-AP» device - the device for maintenance of high factor of health. Its action is based on emission of electromagnetic waves, which resonate with (self) oscillations of parasites of different types - from viruses to helminthes. Thus, DeVita-AP selectively and deliberately influences causal organisms and destroys them, which allows the human body to function in the natural conditions not weighted down by external factors.

The outstanding effect of «DeVita-AP» use is achieved owing to high penetrating power of electromagnetic field into a human body. Frequency of oscillations is set with high precision that allows for influencing a specific microorganism deliberately, without causing any damage to surrounding tissues and organs.

Specially developed care programs allow for maintaining a high health factor in conditions of the most diverse external and internal factors' influence. Each developed course is aimed at influencing a certain microorganism.

The most important advantage of «DeVita-AP» is availability of automatic, specially developed programs for maintenance of high health factor in conditions of the most diverse external and internal factors' influence. Moreover, the owner of the device has a possibility of device programming using any available preventive programs with «DeInfo» utilization. That enables each person to select an individual course of preventive care procedures.

## List of programs

- |                      |                           |
|----------------------|---------------------------|
| 1. Immune resistance | 11. Anti-allergy          |
| 2. Clean eyes K      | 12. No cold               |
| 3. Clean eyes A      | 13. Cocci stop            |
| 4. Eyelids           | 14. Clean eyes            |
| 5. Cornea            | 15. Free breath           |
| 6. Anti-flu S        | 16. Helping a gourmand    |
| 7. Anti-flu B        | 17. Antiseptic            |
| 8. Anti-flu A        | 18. Clean feet            |
| 9. Anti-flu B        | 19. The head free of pain |
| 10. Detoxication     | 20. No worms              |

- |                                |                       |
|--------------------------------|-----------------------|
| 21. Intimate health            | 26. No herpes         |
| 22. Clean lymph                | 27. No protists       |
| 23. Clean liver                | 28. Candida stop      |
| 24. Clean kidneys              | 29. Colon bacillus    |
| 25. Power of liver and kidneys | 30. Anti-inflammation |

**Attention!** The device has contraindication to its use by women in the first trimester of pregnancy and by people, who undergone an organ transplantation. In addition, it is not recommended to use the device within the first two months after suffering a myocardial infarction. In case of suffering from serious heart disease, the device must be positioned no closer than 0.5 m from the body.

### 1. Immune resistance (Imm. resistance)

Immune resistance is an ability of the body to resist the influence of pathogenic viruses, bacteria, and other harmful factors of environment. The stronger the immune resistance is, the less the person is exposed to risk of disease. Strengthening of immune powers is one of the top priorities of Wellness trend. Kicking the bad habits, using balanced food, leading a healthy lifestyle contributes to the strengthening. The Immune resistance program stimulates immune forces of the body and provides reliable protection against external influences.

**Procedure:** the care program needs to be applied daily in the morning, right after waking. It is not recommended to apply the program before going to bed as it can cause sleeplessness.

### 2. Clean eyes K (Clean eyes For Kids)

Children, when playing in a sandbox, bathing in water, communicating with animals may rub their eyes with dirty hands. It sometimes leads to causative agents' of a bacterial conjunctivitis getting into eyes. The Clean eyes D program is aimed at prevention, inflammation removal, killing the pain and smarting in children eyes, blood circulation restoring.

**Procedure:** the program is applied within a 10-day course of procedures, 2-3 times a day. It is recommended to use this program in combination with the Detoxication program once in 3 days.

### 3. Clean eyes A (Clean eyes Anti-Allergy)

During the spring-and-summer period, some people cannot go outdoors because of dust, pollen emerging in the air, which causes blear-eyedness, an intolerable itch, eyes reddening and puffiness. All these symptoms are manifestation of the allergic conjunctivitis occurring also when consuming certain food, medicines or as a reaction to a contact with animal hair. The Clean eyes A program will help to endure all these symptoms easier.

**Procedure:** it is recommended to apply the program 2-3 times a day. It is recommended to use this program in combination with the Detoxication program once in 3 days.

#### 4. Eyelids (Eyelids)

Eyelids – eyes' protection. Moreover, sight condition in many respects depends on eyelids' health. The swelling of eyelids, burning, formation of ulcerations on the brink of eyelids, and eyelashes growth disorder cause a serious discomfort. The outer appearance of the eyes is unappealing and a person starts to suffer from photophobia. The use of Eyelids program brings significant relief in combination with drug treatment, reduces or completely removes the symptoms of the disease.

**Procedure:** the program is applied 2 times a day for 20 days. It is recommended to use this program in combination with the Detoxication program once in 3 days.

#### 5. Cornea (Cornea)

Cornea - an external cover of an eye, which protects it from damages and dust. Keeping the cornea healthy means saving the whole eye. The possible causes of cornea damage are various and can lead to opacification, eyes reddening, puffiness, blur-eyedness. The Cornea program helps getting rid of these phenomena.

**Procedure:** the program is applied 2 times a day for 20 days. It is recommended to use this program in combination with the Detoxication program once in 3 days.

#### 6. Anti-flu S (Anti-flu Swine)

The outbreak of so-called «pig influenza» virus occurred in 2009. It is a virus of the H1N1 subtype, which is transferred from a person to person and which causes diseases with symptoms, typical for flu. During the epidemic period, it is necessary to take preventive measures, the range of which needs to include an application of the Anti-flu C.

**Procedure:** for the purpose of disease prevention during the epidemy, the program is to be applied once in 2 days. If necessary - daily. It is recommended to use this program in combination with the Detoxication program once in 3 days.

#### 7. Anti-flu B (Anti-flu B)

The flu B virus propagates only in human environment. The most widespread way of transmission of the infection is an airborne way. Timely prevention makes it possible to reduce the risk of disease emergence.

**Procedure:** for the purpose of disease prevention during the epidemy, the program is to be applied once in 2 days. If necessary - daily. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 8. Anti-flu A (Anti-flu A)

The virus of A flu is dangerous not only for people, but also for animals. It causes heavy epidemics and covers large areas. Therefore it is necessary to conduct prevention procedures in due time. Use the Anti-flu A program for this purpose. Also, use programs for lymph, liver and kidneys cleaning. Observance of personal hygiene rules of is mandatory.

**Procedure:** for the purpose of disease prevention during the epidemic, the program is to be applied once in 2 days. If necessary - daily. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 9. Anti-flu B (Anti-flu Basic)

Flu is a heavy virus infection. However, it is not the virus itself, but the complications that are dangerous. In conditions of a weakened immune resistance, the body becomes exposed to «attacks» of other infections. To avoid unpleasant consequences during a mass epidemy of the disease, use the Anti-flu B program. For the purpose of bacterial infection prevention, it is also necessary to apply the Cocci stop program and programs for lymph, liver and kidneys cleaning.

**Procedure:** for the purpose of disease prevention during the epidemic, the program is to be applied once in 2 days. If necessary - daily. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 10. Detoxication (Detoxication)

The program assists a liver in intensive struggle with parasites and infections, thus reducing the stress of this important organ and accelerating the excretion of harmful substances.

**Procedure:** the program is applied once in 3 days. When using the program, it is necessary to drink pure drinking water in the amounts of not less than 1-1.5 l a day.

## 11. Anti-allergy (Anti-allergy)

It is impossible to get rid of allergy entirely and simply to reduce its symptoms expression is a difficult task. The Antiallergy program reduces the severity of allergic manifestations - itch, cold, blear-eyedness, reduces red spots on skin. It is very important to increase immune resistance at the same time.

**Procedure:** the program is applied daily in the morning and in the evening during 10 sessions. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 12. No cold (No cold)

The majority of acute respiratory diseases is propagated in airborne way. On a workplace, in public transport, at the cinema – as soon as someone sneezed, all the people around happen in a risk zone at once! Everyone, who stayed in one room with an infected person or got a heavy hypothermia, has to take prevention measures - to use the No cold program. The program will prevent the viruses' propagation and protect the human body.

**Procedure:** for the purpose of disease prevention, the program is to be applied twice a day. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 13. Cocci stop (Cocci stop)

Staphylococci and streptococci are constantly present on the skin and mucous membranes of throat and nose. Often during the trips, people do not get enough sleep, get hypothermia, and get physically tired. Immune resistance weakens. In this case, the cocci already start to have harmful effects on human body. One gets a sore throat, heat, fever. The Cocci stop program reduces negative influence of cocci, enables the human body to restore its powers and immune resistance.

**Procedure:** the program should be applied daily, until a complete recovery of health. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 14. Clean eyes (Clean eyes)

Pain, puffiness, eyes reddening, blear-eyedness may be the expressions of eyes inflammation. Reasons are different. Use the Clean eyes program for eyes condition improvement.

**Procedure:** the program is applied daily as necessary. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 15. Free breath (Free breath)

Free breath is the main constituent of good health. In case of breath disorder the body does not obtain sufficient amount of oxygen, thus the energy, which is so essential for normal body functioning and all the vital processes, is not produced. The Free breath program addresses the issue of nasal blockage caused by various reasons.

**Procedure:** the program is applied once in 2 days. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 16. Helping a gourmand (Help a gourmand)

One often wants to taste various cuisine dishes, which are offered so widely by cafe and restaurants. Nevertheless, the body not always manages to digest all the food, which leads to intestinal upsets, bellyache, nausea, heat and other unpleasant symptoms emergence. The Helping a gourmand program normalizes intestines and stomach functioning, helps the body to excrete dangerous substances.

**Procedure:** it is recommended to apply the program 2-3 times a day until complete recovery of health. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 17. Antiseptic (Antiseptic)

Intense human life is not devoid of small grazes, cuts. Often they are simply not paid attention to. However, in case of solution of skin integuments continuity, infections and various microorganisms, which cause inflammatory process, may penetrate into body. In this case, it is necessary to cleanse them and to apply any anti-septic to the wound. In addition, the Antiseptic program will assist the excretion of harmful substances and quick recovery.

**Procedure:** the program is applied in case of any inflammatory process occurrence. Frequency of application depends on extent of inflammation manifestation (2-3 times a day). After the program completion, it is necessary to use the programs for lymph, liver and kidneys cleaning.

## 18. Clean feet (Clean feet)

It is possible to catch a foot fungus in public places (in pools, at gyms, saunas etc.). The Clean feet program makes a body more resistant against the influence of fungi, reduces their development. Use this program as a prevention of fungoid diseases and you will always have clean and beautiful feet.

**Procedure:** the program should be applied 2 times a day until the symptoms disappear. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 19. The head free of pain (Head w/o pain)

The headache can occur at any time. No need to bear it! When using The head free of pain program a person feels the pain gradually disappearing.

**Procedure:** apply as necessary at any time.

## 20. No worms (No worms)

Most frequently, the infestation with various intestinal infections occurs by way of eating unwashed fruit, meat and fish, which were not thermally conditioned appropriately, contact with dirty hands. Helminthes (worms) can get penetrate into

human organism. When living in a human body the worms eat human tissues, discharge toxic substances.

When undergoing the procedures of No worms program human body is cleaned of these organisms. The program does not exercise any toxic effect on a human body.

**Procedure:** for the preventive purposes, the program is to be applied once in 3 days. For not less than 10 sessions in total. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 21. Intimate health (Intimate health)

Means of protection from sexually transmitted infections are not always able to ensure full safety. To protect oneself, it is possible to use the Intimate health program, which is focused on prevention of being infected with wide range of causative agents of venereal infections.

**Procedure:** it is recommended to apply the program daily during 10 - 20-day course of procedures. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 22. Clean lymph (Clean lymph)

One of the main programs for body cleaning is Clean lymph. Lymphatic system not only performs the function of neutralization of harmful and toxic substances, but also cleans the organs and tissues of them. In cases of suffering from various diseases, the content of harmful substances in human body increases many times. Lymphatic system is not succeeding in excreting them, and recurrent intoxication may occur. The Clean lymph program is designed to accelerate the lymph cleaning, while not stressing the lymphatic system itself. It is recommended to apply the program one week prior to a vacation trip, which will help to improve immune resistance and to prepare the body for acclimatization.

**Procedure:** it is recommended to apply the program once in 3 days. When using the Clean lymph program, it is necessary to drink pure drinking water in the amounts of not less than 1-1.5 l a day.

## 23. Clean liver (Clean liver)

The liver is a main organ, which neutralizes toxic and other harmful substances, which form in human body. Thus, the liver is heavily stressed, its working capacity decreases. It is necessary to use the program as a preventive measure, especially after undergoing the programs of parasites and infections counteracting.

**Procedure:** the program is applied once in 3 days. When using the Clean liver program, it is necessary to drink pure drinking water in the amounts of not less than 1-1.5 l a day.

## 24. Clean kidneys (Clean kidneys)

Kidneys clean the body from the water and the substances dissolved in it, which were discharged by organs and tissues, but are not needed any more. It is very important to assist in cleaning of kidneys, especially during infections outbreak, to enable them excrete the useless and harmful substances from the body. This specific task is performed by the Clean kidneys program.

**Procedure:** the program is applied once in 3 days. When using the Clean kidneys program, it is necessary to drink pure drinking water in the amounts of not less than 1-1.5 l a day.

## 25. Power of liver and kidneys (Liver & kidneys)

The program should be used when feeling sick (when the liver and kidneys function in a stressed mode to neutralize and excrete harmful and toxic substances), and after undergoing the programs for liver and kidneys cleaning. The program is aimed at maintaining the performance capability of these organs and restoration of their abilities.

**Procedure:** the program is applied once in 3 days. When using the Power of liver and kidneys program, it is necessary to drink pure drinking water in the amounts of not less than 1-1.5 l a day.

## 26. No herpes (No herpes)

"Cold sores on lips" is a most widespread form of herpes. It can also manifest itself on genitals, neck, head and hands, in the form of small ulcerations and bladders. Herpes infection can penetrate into nervous system causing its various disorders. When inpouring into human body, herpes penetrates into nervous cells and remains there, and then manifests itself only in case the immune resistance is lowered. The carrier of the virus can infect other people via kisses, sexual contact and even through exposed regions of skin integument.

**Procedure:** it is recommended to use the No herpes program 2-3 times a day for one week, and then once a day for 3 weeks. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 27. No protists (No protists)

The program prevents the negative impact of a large number of monocellular parasites on a human body. They inpour into body through the contact with dirty hands, consuming unwashed food and unboiled water. The program affects the protists, destroying them. Strict observance of personal hygiene rules and regular use of No protists program as preventive measures protect the person from these microorganisms.

**Procedure:** it is recommended to use the program 1-2 times a day for 10-14 days. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 28. Candida stop (Candida stop)

The program protects a person from infecting with such widespread fungus, as Candida. This microorganism is a part of normal microflora of intestines, vagina, oral cavity. However, when the immune resistance is lowered, the Candida starts to breed intensely and harm the body. Candidiasis prevention by means of Candida stop program normalizes acid-base balance, which restrains fungus propagation.

**Procedure:** the program should be applied 1-3 times a day for 2-4 weeks. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 29. Colon bacillus (Colon bacillus)

The program is capable of suppressing negative effect of Colon bacillus and increase body resistance to it. Pathogenic Colon bacilli, when getting into a human body with low quality and badly cooked food or water, cause severe intoxication. Use of the Colon bacillus program helps protecting the body from harmful microorganisms, while not damaging the intestines microflora.

**Procedure:** it is recommended to apply the program daily for 10-20 days. It is recommended to use this program in combination with the Detoxication program once in 3 days.

## 30. Anti-inflammation (Anti-inflammat.)

Inflammation can emerge at various trauma conditions or upon infection's penetration into the body. The Anti-inflammation program helps in restoring blood circulation, in removing various toxic substances, which are forming during an inflammation. Thus, tissues oxygen supply and metabolic processes improve, which aids quicker recovery.

**Procedure:** the program should be applied under various inflammatory conditions for the purpose of their reduction. Frequency of application is 1-3 times a day. For the purpose of detoxication after completing this program, it is necessary to apply programs of lymph, liver and kidneys cleaning.